



*Sticker for illustration only, does not include AED.



ADDITIONAL INFORMATION

Weight	5 kg
Dimensions	38 × 44 × 19 cm
Pricing Unit	Kit
Size	44 x 38 x 15.4cm
Min Order Qty	1
ARTG ID	406543
Medical Device Class	Class 2a
Tax Schedule	Taxable
Item Length	44
Item Width	38
Item Height	15.4
Item Measurement Unit	cm
Location	AU-Cameron Dr
Is Serialized Item	No
Kits - People (High Risk)	10
Kits - People (Low Risk)	25
Kits - Treatments	Burns, CPR, Dressings and Bandages, Eye Injury, Major Wounds, Minor Wounds, Snake Bite
Kits - Industry	Construction / Outdoor / Remote Location
Kits - Case	Soft Pack
Is Dangerous Good	No

PRODUCT DESCRIPTION

The AeroKit Modulator Extreme waterproof tough first aid and trauma kit is designed to give rapid, organised support when the unexpected happens. Built around a heavy duty, high visibility case, it helps protect critical first aid supplies from dust, water, and rough handling so they are ready when needed. Inside, colour coded injury specific modules group items for cuts and grazes, burns, eye injuries, dressings and bandages, snake bite, trauma and advanced bleed control. This structure lets users move straight to the right module with little more than a glance, which is especially valuable in high pressure moments at work, at home, or outdoors.

The kit layout supports clear visibility of contents when the lid is open, with modules that can be lifted out and taken directly to the

patient. A dedicated cavity is included for an automated external defibrillator so that an AED such as LifePak CR2 or Heartsine devices can be housed alongside bleed control and CPR supplies. This makes the Modulator Extreme well suited to workplaces, vehicles, remote properties and community locations where a combined first aid and defibrillation point is required. The robust case features a carry handle for easy transport around large sites or to incident scenes.

Each module is clearly labelled for fast identification and colour coded to match common injury types, helping less experienced first aiders find the right items quickly. The modular layout also makes it easier to check stock levels and restock specific injury areas without emptying the entire kit. Combined with bleed control components such as a RapidStop tourniquet, trauma dressings and Israeli bandages, the kit supports response to both everyday injuries and more serious events.

Compact external dimensions around 44 cm x 38 cm x 15.4 cm make the kit practical to store in vehicles, workshops, offices and site sheds while still providing generous internal capacity. The hard plastic case is waterproof to an IP67 rating and crush resistant, helping protect the contents from rain, hoses, dust, and rough transport in utes or service vehicles. High visibility yellow colouring assists with quick location in low light or crowded workplaces and also suits home garages and hobby spaces where clear identification of emergency gear is important.

- **Organised response:** colour coded injury specific modules support fast selection of the right first aid items.
- **Rugged protection:** hard plastic IP67 rated case resists water, dust and crushing for dependable performance.
- **Combined first aid and trauma:** includes space for an AED alongside modules for cuts, burns, eye injuries, trauma and snake bite.
- **Easy to carry:** integrated handle and balanced case design support quick grab-and-go use around workplaces, homes and remote sites.
- **Ideal for work, home and travel:** suitable for vehicles, workshops, farms, construction sites and household emergency stations.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.