

AEROHAZARD NEEDLE DISPOSAL UNIT 90ML

ND90



ADDITIONAL INFORMATION

Weight 0.025 kg

Dimensions $2.5 \times 6 \times 9 \text{ cm}$

Carton Qty 200

Pricing Unit Unit/1

Size 90ml

Min Order Qty 1

ARTG ID N/A

GMDN Code N/A

UNSPSC 42192802

Tax Schedule Taxable

Is Medical Device No

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

Aero offers a wide range of waste disposal containers designed to meet user demands for strength, simplicity and value. This range offers outstanding protection against penetrating sharps and leaking fluids due to its rigid polypropylene construction and wall thickness.

Needle Disposal Units are a durable and convenient locking small sharps disposal safe for needles, syringes and scalpel blades.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

AeroHazard Needle Disposal Unit



AEROHAZARD NEEDLE DISPOSAL UNIT 90ML

ND90