



ADDITIONAL INFORMATION

Weight	0.025 kg
Dimensions	2.5 × 6 × 9 cm
Carton Qty	200
Pricing Unit	Unit/1
Size	90ml
Min Order Qty	1
ARTG ID	N/A
GMDN Code	N/A
UNSPSC Code	42192802
Tax Schedule	Taxable
Is Medical Device	No
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

Aero offers a wide range of waste disposal containers designed to meet user demands for strength, simplicity and value. This range offers outstanding protection against penetrating sharps and leaking fluids due to its rigid polypropylene construction and wall thickness.

Needle Disposal Units are a durable and convenient locking small sharps disposal safe for needles, syringes and scalpel blades.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AEROHAZARD NEEDLE DISPOSAL UNIT 90ML

ND90