



#### **ADDITIONAL INFORMATION**

Weight 0.14 kg

**Dimensions**  $19.5 \times 4 \times 8.5 \text{ cm}$ 

Carton Qty 24

Pricing Unit Unit/1

Size Adult

Min Order Qty 1

**ARTG ID** 145376

GMDN Code 18009

Medical Device Class 2a

**UNSPSC** 42142620

**Tax Schedule** Taxable

**Is Medical Device** Yes

**Is Serialized Item** No

Is Dangerous Good No

#### **PRODUCT DESCRIPTION**

#### NIO Adult = New Intraosseous Device

## NO DRILLS, BATTERIES, OR EXTRA PARTS. MEET THE NEW IO.

NIO offers unparalleled performance, ergonomics and usability for rapid and safe intraosseous access to multiple locations.

NIO is a simple, automatic IntraOsseous (IO) device ready to use for safe, quick and easy vascular access.

NIO is used in emergency situations when vascular access is needed for fluid and drug administration. NIO is small, light, and easy to use in even the most demanding situations.

NIO is designed to be easily injected into either the tibia or the humerus in all clinical settings.

Features:



## NIO ADULT WITH FIXATION DRESSING

NIO-A

- Sterile
- o Ready To Use
- User-friendly
- No exposed needle
- o Pocket-sized & lightweight ~100g
- Needle Gauge 15G
- No external power source or battery
- Vascular access in 10 seconds
- Disposable and single use
- No maintenance required
- Safety features
- o 5-year shelf life

### **CLINICALLY PROVEN SUCCESS**

#### SAFE ON THE BONE.

The NIO has proven to be safe on the bone and resulted in zero incidents of bone fracture as observed by post-insertion x-rays.

#### **EASILY TOLERATES COMMON USER ERRORS.**

A recent study evaluated common user errors related to needle placement and found that even when placed up to 5cm from the correct insertion site, the NIO maintained an 85% success rate.

#### FIRST ATTEMPT SUCCESS RATES AS HIGH AS 98%.

Initial studies demonstrated an overall average success rate of 91.3% for the proximal tibia and 93.1% for the humeral head. The New York Institute of Osteopathic Medicine conducted an additional study which established a proximal tibia success rate of over 98%.

"...the first attempt success rate for inserting an NIO in the tibia was 98%, Humerus 91%..."

NEW YORK INSTITUTE OF TECHNOLOGY COLLEGE OF OSTEOPATHIC MEDICINE



# NIO ADULT WITH FIXATION DRESSING

NIO-A

## **Product Specifications**

139mm (5.47in) L

Dimensions 64.5mm (2.6in) W

43.7mm (1.72in) H

Weight 94.5g

Packaging Packed individually in a peel-back blister pack

Needle Gauge 15G

Needle Effective Length . 42mm (1.65in) Needle Penetration Depth 25mm (0.98in)

1. Rotating safety latch prevents accidental activation

Safety Mechanism

2. Simultaneously compressing the device against the skin while squeezing the trigger wings

activates device

Shelf Life\* 5 years

Specifications are subject to change without notice.

\*From date of manufacture

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

NIO Instraosseous device