

# NIO SIMULATION IO TRAINING DEVICE PAED (AGE 3-12)

NIO-T-P



#### **ADDITIONAL INFORMATION**

Weight 0.12 kg

**Dimensions**  $18 \times 4.9 \times 9 \text{ cm}$ 

Carton Qty 12

**Pricing Unit** Each

Min Order Qty 1

Tax Schedule Taxable

Is Serialized Item No
Is Dangerous Good No

### PRODUCT DESCRIPTION

### NIO Adult = New Intraosseous Device

NIO offers unparallel performance, ergonomics and usability for rapid and safe intraosseous access to multiple locations.

NIO is a simple, automatic IntraOsseous (IO) device ready to use for safe, quick and easy vascular access.

NIO is used in emergency situations, when vascular access is needed for fluid and drug administration. NIO is small, light, and easy to use in even the most demanding situations.

NIO is designed to be easily injected into either the tibia or the humerus in all clinical settings.

#### Features:

- Sterile
- · Ready To Use
- User friendly
- No exposed needle
- Pocket sized & lightweight ~100g



# NIO SIMULATION IO TRAINING DEVICE PAED (AGE 3-12)

NIO-T-P

- No external power source or battery
- Vascular access in 10 seconds
- Disposable and single use
- No maintenance required
- Safety features
- 5 year shelf life

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.