



## ADDITIONAL INFORMATION

<b>Weight</b>	0.08 kg
<b>Dimensions</b>	17.4 × 4 × 7.5 cm
<b>Pricing Unit</b>	Each
<b>Min Order Qty</b>	1
<b>Tax Schedule</b>	Taxable
<b>Location</b>	AU-Warehouse
<b>Is Serialized Item</b>	No

## PRODUCT DESCRIPTION

The NIO Infant Simulation device provides simulated and realistic training for safe use of NIO Infant devices.

The NIO Infant Simulation includes a Fixation dressing, which secures the stabilizer base to the skin and prevents dislodgment once intraosseous access has been achieved.

The NIO Infant is a ready to use device and requires no batteries, no assembly, and no extra parts.

Can be used on human skin without penetration, suitable for cost effective training of intraosseous device application.