



ADDITIONAL INFORMATION

Weight	0.7 kg
Dimensions	24 × 22 × 9 cm
Carton Qty	1
Pricing Unit	Each
Min Order Qty	1
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The NIO Infant™ Simulation offers an unrivalled IO simulation experience. Its training leg includes two interchangeable bone models for varying bone density, allowing users to learn and experience first-hand an IO device.

FEATURES:

- NIO-I SIM
- Infant Training Leg
- (x3) NIO Infant Fixation demo dressing (non-adhesive)
- Bone model white - simulates young infant bone
- Bone model pink - simulates older infant/toddler bone
- The bone model within the training leg can be rotated to allow for additional training attempts. The bone model is suitable for fluid administration