

# **NIO INFANT SIM (TRAINER) KIT**

NIO-T-IK



### **ADDITIONAL INFORMATION**

Weight 0.7 kg

**Dimensions**  $24 \times 22 \times 9 \text{ cm}$ 

Carton Qty 1

Pricing Unit Each

Min Order Qty 1

Tax Schedule Taxable

Is Serialized Item No
Is Dangerous Good No

## PRODUCT DESCRIPTION

The NIO Infant™ Simulation offers an unrivalled IO simulation experience. Its training leg includes two interchangeable bone models for varying bone density, allowing users to learn and experience first-hand an IO device.

#### **FEATURES:**

- NIO-I SIM
- Infant Training Leg
- (x3) NIO Infant Fixation demo dressing (non-adhesive)
- Bone model white simulates young infant bone
- Bone model pink simulates older infant/toddler bone
- The bone model within the training leg can be rotated to allow for additional training attempts. The bone model is suitable for fluid administration

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



# NIO INFANT SIM (TRAINER) KIT

NIO-T-IK