

# **NIO PLUS ADULT WITH NIO FIXATION DRESSING**

NIO-PL-A



## ADDITIONAL INFORMATION

Waight	0.16 kg
Weight	0.16 kg
Dimensions	16 × 7 × 4 cm
Carton Qty	12
Pricing Unit	Each
Min Order Qty	1
UNSPSC	42142620
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

#### **PRODUCT DESCRIPTION**

The NIO+ Adult provides safe, rapid intraosseous access in the most extreme environments under extreme conditions. The NIO+ Adult is a reliable, automatic, single-use intraosseous device, robustly manufactured with a reinforced polymer core, which resists elevated thermal temperatures, enabling the device to operate in extreme environments. Light weight, at only 140grams the device can be stored within a small medical pouch, making it a practical solution for emergency situations.

- Compact and durable packaging
- Proven durability in extreme temperatures
- 5-year room temperature and 3-year extreme temperature shelf life
- Vascular access in 10 seconds
- Needle Gauge 15G
- Most cost-effective IO solution



**NIO-PL-A** 

- Lightweight approx. 140g (packaged)
- Dual safety features
- Semi-automatic device
- No exposed needle
- Single-use device
- Fully assembled

# **Product Specifications**

Dimensions	139 mm (5.47") L 57 mm (2.24") W 27 .33 mm (1.06") H	
Weight	90g (140g packaged)	
Packaging	Packed individually	
Needle Gauge	15G	
Needle Effective Length .	42mm (1.65in)	
Needle Penetration Depth 25mm (0.98in)		
Safety Mechanism	<ol> <li>The cap must be rotated prior to deployment.</li> <li>Device activation requires the user simultaneously press the device to the insertion site and pull its trigger.</li> </ol>	
Shelf Life*	3 years up to 50°C or 5 years at room temperature	

## \*From date of manufacture

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

NIO Instraosseous Device