



ADDITIONAL INFORMATION

Weight	0.16 kg
Dimensions	16 × 7 × 4 cm
Carton Qty	12
Pricing Unit	Each
Min Order Qty	1
UNSPSC Code	42142620
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

The NIO+ Adult provides safe, rapid intraosseous access in the most extreme environments under extreme conditions. The NIO+ Adult is a reliable, automatic, single-use intraosseous device, robustly manufactured with a reinforced polymer core, which resists elevated thermal temperatures, enabling the device to operate in extreme environments. Light weight, at only 140grams the device can be stored within a small medical pouch, making it a practical solution for emergency situations.

- Compact and durable packaging
- Proven durability in extreme temperatures
- 5-year room temperature and 3-year extreme temperature shelf life
- Vascular access in 10 seconds
- Needle Gauge 15G
- Most cost-effective IO solution

- Lightweight approx. 140g (packaged)
- Dual safety features
- Semi-automatic device
- No exposed needle
- Single-use device
- Fully assembled

Product Specifications

Dimensions	139 mm (5.47") L 57 mm (2.24") W 27 .33 mm (1.06") H
Weight	90g (140g packaged)
Packaging	Packed individually
Needle Gauge	15G
Needle Effective Length	42mm (1.65in)
Needle Penetration Depth	25mm (0.98in)
Safety Mechanism	1.The cap must be rotated prior to deployment. 2. Device activation requires the user simultaneously press the device to the insertion site and pull its trigger.
Shelf Life*	3 years up to 50°C or 5 years at room temperature

*From date of manufacture

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

NIO Intraosseous Device