



ADDITIONAL INFORMATION

Weight	0.36 kg
Dimensions	21.5 × 3.5 × 3.5 cm
Carton Qty	1
Pricing Unit	Each
Min Order Qty	1
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

- Simulates older infant (pink) and younger infant (white) bone
- Rotatable for multiple training attempts
- Bone density allows for realistic feel for users
- Suitable for fluid administration

NIO SIM Infant bone models: Bone Model white (1) / Bone model pink (1)