



ADDITIONAL INFORMATION

Weight 0.08333 kg

Carton Qty 5

Pricing Unit Unit/1

Size One Size

Min Order Qty 1

Tax Schedule Taxable

Is Serialized Item No
Is Dangerous Good No

PRODUCT DESCRIPTION

NIO Adult = New Intraosseous Device

NIO offers unparallel performance, ergonomics and usability for rapid and safe intraosseous access to multiple locations.

NIO is a simple, automatic IntraOsseous (IO) device ready to use for safe, quick and easy vascular access.

NIO is used in emergency situations, when vascular access is needed for fluid and drug administration. NIO is small, light, and easy to use in even the most demanding situations.

NIO is designed to be easily injected into either the tibia or the humerus in all clinical settings. Features:

- Sterile
- Ready To Use
- o User friendly



NIO TRAINER RELOADER

NIO-T-R

∘ No expo	osed needle		
○ Pocket s	sized & lightweight ~100g		
No exte	ernal power source or battery		
o Vasculai	or access in 10 seconds		
∘ Disposa	able and single use		
○ No mair	ntenance required		
○ Safety fe	eatures		
∘ 5 year s	shelf life		



NIO TRAINER RELOADER

NIO-T-R

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.