



ADDITIONAL INFORMATION

Weight	0.54 kg
Dimensions	15.8 × 13.7 × 13.8 cm
Carton Qty	8
Pricing Unit	Box/50
Size	10cm x 20cm
Min Order Qty	1
ARTG ID	159437
GMDN Code	34654
Medical Device Class	Class 1s
UNSPSC Code	42311545
Tax Schedule	Taxable
Item Length	20
Item Width	10
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No
Shelf Life (Months)	60
Warranty Period (Months)	0



PRODUCT DESCRIPTION

AeroPad™ Non-Adherent Dressing is a soft, non-stick wound pad designed to protect delicate skin while offering reliable absorbency. A



AEROPAD NON-ADHERENT DRESSING 10 X 20CM BOX/50

APD102S

dual-sided non-adherent contact layer helps minimise trauma during dressing changes, while the cushioned inner core draws fluid away from the site to support comfort. The dressing is sterile and latex free, making it suitable for use across workplaces, schools, clinics and home first aid kits.

The pad sits comfortably beneath securement such as film dressings, paper tapes or conforming bandages, allowing teams to build flexible wound-cover solutions for cuts, grazes and minor burns. Its low-lint, non-woven construction conforms to body contours and is easy to position for effective coverage on hands, joints and other high-movement areas.

Engineered for dependable performance, AeroPad™ helps standardise wound care across busy environments. The combination of a gentle non-adherent surface and an efficient absorbent layer supports practical changeovers and user comfort, making it a smart choice for day-to-day first aid readiness and clinical triage.

- **Non-adherent contact layer:** helps reduce pain and trauma on removal.
- **Absorbent inner core:** manages exudate for drier, more comfortable wear.
- **Sterile and latex free:** suitable for sensitive users and diverse settings.
- **Soft, conformable pad:** low lint design sits neatly under tape or film.

The **Size 10 x 20 cm** pad provides broad coverage for larger areas or post-procedure sites while maintaining a gentle, non-stick contact layer and soft absorbency for comfortable wear.

Supplied in a **box of 50** individually pouched, sterile pads to streamline replenishment for clinics, workplaces and first aid stations while supporting efficient everyday care.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Low Adherent Dressing Pad Box/50