

OXYGEN THERAPY MASK WITH 2M TUBING -**ADULT**

OX003



ADDITIONAL INFORMATION

Weight 0.06 kg

Dimensions 12.9 × 13.9 × 7.6 cm

Carton Qty 40

Pricing Unit Each

Size Adult

Min Order Qty 1

ARTG ID 224158

GMDN Code 42424

Medical Device Class Class 1

UNSPSC 42271708

Tax Schedule Non-Taxable

Is Medical Device Yes Is Serialized Item No

No

Is Dangerous Good

PRODUCT DESCRIPTION

The AeroRescue Oxygen Therapy Mask supports reliable delivery of supplemental oxygen in first response, workplace first aid and clinical settings. The transparent facepiece allows easy visual assessment, while the contoured shape and soft edges help achieve a comfortable seal over the nose and mouth. An adjustable elastic strap assists with quick placement and secure positioning during treatment, even when the patient or responder is in motion.

The mask interfaces with standard oxygen equipment used in emergency care. Its clear construction aids communication and observation, and the low-profile design helps maintain access to the airway and surrounding facial landmarks for monitoring. Built for compatibility with common oxygen regulators and tubing, it integrates smoothly into resuscitation workflows alongside bag valve masks and basic life support equipment.

Designed to assist responders in delivering air/oxygen gases to a patient, this mask is a practical choice for first aid rooms, ambulance kits and training environments. It provides dependable performance during oxygen therapy and resuscitation support without adding complexity to time-critical care.

- Secure, comfortable fit: contoured mask with adjustable strap for quick, stable placement.
- Clear visibility: transparent body supports patient assessment and communication.
- Equipment ready: connects with common oxygen delivery setups used in first response.
- Workflow friendly: low-profile design supports efficient airway management.



OXYGEN THERAPY MASK WITH 2M TUBING - ADULT

OX003

This prod	duct is not	a substitute	for direct me	dical attent	ion from a	healthcare	practitioner.	Always read	the label	and follow	the t
direction	is for use.										

Oxygen Therapy Mask & Tubing 2m