



ADDITIONAL INFORMATION

Weight	0.011 kg
Dimensions	10 × 0.3 × 14 cm
Pricing Unit	Qc-emsr each
Size	7.5cm x 1.2M
Min Order Qty	1
ARTG ID	314283
GMDN Code	47917
Medical Device Class	Class 2b
UNSPSC Code	51131800
Tax Schedule	Taxable
Item Length	120
Item Width	7.5
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

QuikClot EMS Rolled Gauze is 3 inches x 4 feet of soft, white, nonwoven, hydrophilic gauze impregnated with kaolin, an inorganic mineral that is both safe and effective in accelerating the body's natural clotting cascade without any exothermic reactions or use of animal or human proteins.

Each QuikClot® EMS Rolled Gauze, is packaged in an easy-tear pouch.

QUIKCLOT® Combat Gauze is the only product carried by all branches of the US Military to control life-threatening hemorrhage.

QuikClot EMS 4x4 Dressing features:

- Effective, Simple, and Safe
- Easy to use and remove
- Flexible and pliable
- 2 Year Expiry from Date of Manufacture

The QuikClot range of hemostatic products have been chosen by EMS, Law Enforcement, Military and Healthcare professionals to help



QUIKCLOT HAEMOSTATIC EMS ROLL 7.5CM X 1.2M

QC-EMSR

save lives everyday.



QUIKCLOT HAEMOSTATIC EMS ROLL 7.5CM X 1.2M

QC-EMSR

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

QuikClot EMS Roll