

ADDITIONAL INFORMATION

Weight	2.76 kg
Dimensions	60.5 × 27.7 × 18 cm
Carton Qty	1
Pricing Unit	Qty 1
Min Order Qty	1
ARTG ID	N/A
GMDN Code	N/A
UNSPSC Code	42301502
Tax Schedule	Taxable
Is Medical Device	No
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

For years, instructors of infant CPR have not been able to know for sure what the participants learned. Until now. Little Baby QCPR is a realistic and affordable BLS manikin for pediatrics. By utilizing gamified learning and objective feedback, Little Baby QCPR improves training quality, learner engagement, and classroom efficiency.

EASY TO USE. QUALITY ASSURED

Set up a class in no-time

Use the one-click connection to set up training. The QCPR app onboarding guides you through the setup of your classroom.

Engage your learners

Increase engagement and confidence in both learner and instructor with a combination of realistic feel and intuitive feedback.

Guideline compliant

With objective feedback on compressions and ventilations, Little Baby QCPR is future-proof. It teaches all the parameters of high-quality CPR as defined by the AHA.

INSTALL CONFIDENCE IN YOUR LEARNERS

Foreign body removal

Choking is one of the biggest reasons for cardiac arrest in infants. With Little Baby QCPR you allow learners to practice the not so gentle touch to see what it takes to remove airway obstructions.

Correct hand placement

Spot incorrect hand placement with built-in sensors and QCPR app to instantly adjust learner's performance.

Act when it matters

Practice correct compressions, ventilations, and foreign body removal so that learners act with confidence when it matters the most.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.