



ADDITIONAL INFORMATION

Weight	0.04 kg
Dimensions	16 × 16 × 2 cm
Pricing Unit	Each
Size	7.5cm x 3.65M
Min Order Qty	1
ARTG ID	168804
GMDN Code	47917
Medical Device Class	Class 2b
UNSPSC	51131800
Tax Schedule	Taxable
Item Length	365
Item Width	7.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

QUIKCLOT® Combat Gauze is a soft, white, sterile, nonwoven 7.5 cm x 3.7 m gauze strip made with kaolin, a chemically inactive mineral that does not contain animal or human proteins or botanicals. Each strip of QUIKCLOT® Combat Gauze is individually wrapped in an easy rip, military grade foil pouch. Indicated for temporary external control of traumatic bleeding, QUIKCLOT® Combat Gauze is flexible, pliable and has the ability to contour to all wounds. Recommended as the number one hemostatic agent by the COTCCC (Committee on Tactical Combat Casualty Care Committee).

QUIKCLOT® Combat Gauze is the only product carried by all branches of the US Military to control life-threatening hemorrhage.

- Effective, Simple, and Safe
- Easy to use and remove
- Flexible and pliable

- 7.5 cm x 3.7 m gauze strip

QuikClot is recommended for:

- Gunshot Wounds
- Stabbings
- Mass Casualty Events
- Large Wounds/ Excessive Bleeding

The QuikClot range of hemostatic products have been chosen by EMS, Law Enforcement, Military and Healthcare professionals to help save lives everyday.



QUIKCLOT HAEMOSTATIC COMBAT GAUZE 7.5CM X 3.65M

QC-CG



QUIKCLOT HAEMOSTATIC COMBAT GAUZE 7.5CM X 3.65M

QC-CG

QuikClot Combat Gauze