



## ADDITIONAL INFORMATION

<b>Weight</b>	0.04 kg
<b>Dimensions</b>	16 × 16 × 2 cm
<b>Pricing Unit</b>	Each
<b>Size</b>	7.5cm x 3.65M
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	314283
<b>GMDN Code</b>	47917
<b>Medical Device Class</b>	Class 2b
<b>UNSPSC Code</b>	51131800
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	365
<b>Item Width</b>	7.5
<b>Item Measurement Unit</b>	cm
<b>Is Medical Device</b>	Yes
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

## PRODUCT DESCRIPTION

QUIKCLOT® Combat Gauze is a soft, white, sterile, nonwoven 7.5 cm x 3.7 m gauze strip made with kaolin, a chemically inactive mineral that does not contain animal or human proteins or botanicals. Each strip of QUIKCLOT® Combat Gauze is individually wrapped in an easy rip, military grade foil pouch. Indicated for temporary external control of traumatic bleeding, QUIKCLOT® Combat Gauze is flexible, pliable and has the ability to contour to all wounds. Recommended as the number one hemostatic agent by the COTCCC (Committee on Tactical Combat Casualty Care Committee).

QUIKCLOT® Combat Gauze is the only product carried by all branches of the US Military to control life-threatening hemorrhage.

### QuikClot EMS 4x4 Dressing features:

- Effective, Simple, and Safe
- Easy to use and remove
- Flexible and pliable
- 7.5 cm x 3.7 m gauze strip
- 2 Year Expiry from Date of Manufacture



# QUIKCLOT HAEMOSTATIC COMBAT GAUZE 7.5CM X 3.65M

QC-CG

## RECOMMENDED USE:

### **Smaller Surface Bleeding (Nuisance)**

- Gunshot Wounds
- Stabbings
- Mass Casualty Events
- Large Wounds/ Excessive Bleeding

The QuikClot range of hemostatic products have been chosen by EMS, Law Enforcement, Military and Healthcare professionals to help save lives everyday.



# QUIKCLOT HAEMOSTATIC COMBAT GAUZE 7.5CM X 3.65M

QC-CG

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

QuikClot Combat Gauze