

# QUIKCLOT HAEMOSTATIC EMS DRESSING 10 X 10CM

QC-4X4D



### ADDITIONAL INFORMATION

Weight	0.007 kg
Dimensions	13.3 × 0.1 × 12.5 cm
Pricing Unit	Unit/1
Min Order Qty	1
ARTG ID	314283
GMDN Code	47917
Medical Device Class	Class 2b
UNSPSC	51131800
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

### **PRODUCT DESCRIPTION**

QuikClot EMS 4x4 Dressing is a sterile, soft white 4 in. x 4 in. nonwoven hemostatic dressing with X-ray indicator, and is impregnated with kaolin, an inorganic mineral that activates Factor XII, which in turn accelerates the body's natural clotting ability. This bleeding control solution creates a robust clot to control bleeding fast potentially giving you more time to attend to other tasks as you move patients through the continuum of care. Compared to standard gauze or other available hemostatic devices, QuikClot has superior results that indicate it can reduce blood loss and promote stronger clot formation. These results may help reduce exposure to blood and bloodborne pathogens, and potentially reduce decontamination time for equipment and vehicles.

### **QuikClot EMS 4x4 Dressing features:**

- Intuitive to use, EMS 4x4 Dressing uses the same placement technique as standard gauze
- Contains an X-ray detectable strip for easy identification
- Dressing may be left in place for up to 24 hours, (refer to the IFU for full Instructions for use)
- 2 Year Expiry from Date of Manufacture

#### **RECOMMENDED USE:**

## **Smaller Surface Bleeding (Nuisance)**

- Head Lacerations
- Bleeding AV Fistulas
- Ruptured Varicose Veins



# QUIKCLOT HAEMOSTATIC EMS DRESSING 10 X 10CM

QC-4X4D

- Nosebleeds
- Skin Tears
- Other Smaller Surface Nuisance Bleedings



# QUIKCLOT HAEMOSTATIC EMS DRESSING 10 X 10CM

QC-4X4D

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.