



ADDITIONAL INFORMATION

Weight	0.1448 kg
Dimensions	5 × 19 × 4 cm
Carton Qty	25
Pricing Unit	Each
Min Order Qty	1
ARTG ID	376893
GMDN Code	14072
Medical Device Class	Class 1
UNSPSC Code	42171506
Tax Schedule	Taxable
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

RapidStop Tourniquet is a purpose-built emergency limb tourniquet for major haemorrhage when seconds count. Its intuitive design enables rapid one-hand self-application and even faster assisted use by a responder. Trusted across defence, public safety, industrial safety and outdoor communities, it delivers dependable performance in high-stress conditions where clarity and speed are vital.

The mechanical advantage system combines a wide pull-strap, high-strength webbing and a precision ratcheting buckle to tighten efficiently with gross motor movements. Designed to achieve full occlusion quickly when correctly applied, it is suitable for use on arms and legs and provides a clear, step-wise operation that supports confident action under pressure.

Built for training and field deployment, the durable components resist demanding environments, while the low-profile form factor stages easily in kits, pouches and vehicles. It integrates smoothly into workplace and community bleeding control programs and is a reliable addition to trauma kits for first responders, safety teams and outdoor users.

- **Fast, intuitive tightening:** supports one-hand application when self-aid is required.

- **Mechanical advantage system:** efficient force with minimal effort for reliable control.
- **Robust build:** durable webbing and hardware for demanding environments.
- **Compact, kit-ready form factor:** easy to stage in pouches and first aid kits.
- **Arms and legs:** designed for major limb bleeding control when correctly applied.

The **Military Black** colour option pairs the RapidStop Tourniquet design with a low-visibility finish that blends with tactical uniforms, body armour and dark gear while retaining the same mechanical advantage system and intuitive one-hand operation.

Supplied as a single unit in the **Military Black** finish, this option suits users who want a discreet tourniquet that integrates with tactical rigs, duty belts or dark-coloured first aid kits. Wide strap and ratcheting buckle for swift tightening.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.