



ADDITIONAL INFORMATION

Weight	3.5 kg
Dimensions	40 × 55 × 17 cm
Carton Qty	2
Pricing Unit	Bag/1
Size	54cm x 48cm x 32cm
Min Order Qty	1
UNSPSC Code	42192420
Tax Schedule	Taxable
Item Length	54
Item Width	48
Item Height	32
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

The AeroCase range of softpacks have been specifically designed to meet the needs that exist within the first aid industry. These bags are ideal for portable usage and where space is limited.

ABLR03

The ABLR03 is a durable, red trauma backpack designed for emergency responders. It features reflective tape for visibility in low-light conditions, comfortable carry handles, and multiple pockets for storing essential equipment. The main compartment has modular dividers for organization and the capacity to hold a size C oxygen cylinder.

Perfect for Oxygen Kits and more, this multipurpose bag is an everyday winner.

Please note that this item is not designed for IVR or Inhalo units.

Also available in a smaller trauma case ([ABLR02](#)).



AEROBAG RED TRAUMA FIRST AID BACKPACK 48 X 54 X 32CM

ABLR03

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Red First Aid Backpack Trauma Large