

RESPONDER 4 SERIES PLASTIC WATERPROOF FIRST AID KIT 36 X 8.5 X 25CM

AFAK4W





ADDITIONAL INFORMATION

Weight 1.54 kg

Dimensions $27 \times 36 \times 10 \text{ cm}$

Carton Qty 4

Pricing Unit Unit/1

Size 36cm W x 8.5cm D x 25cm H

Min Order Qty 1

ARTG ID 406543

Tax Schedule Taxable

Item Length 36

Item Width 8.5

Item Height 25

Item Measurement Unit cm

Is Serialized Item No

Kits - People (High Risk) 10

Kits - People (Low Risk) 25

Burns, CPR, Dressings and Bandages, Eye Injury,

Kits - Treatments Major Wounds, Minor Wounds, Snake Bite, Spider

Bite

Construction / Outdoor / Remote Location, Education

Kits - Industry – Grounds Keeping, Education – Office / Classroom /

Hall, Office / Retail

Kits - Case Weatherproof

Is Dangerous Good No

PRODUCT DESCRIPTION

Ensure workplace compliance for first aid with this comprehensive range, containing a wide scope of products to cater for any workplace emergencies. Respond promptly to incidents with a variety of outer cases to suit all types of environments. Internal compartments and shelves keep product organised, and all plastic and metal case options are able to be wall-mounted. The soft pack style is ideal for mobile workers, while the durable plastic cases are recommended for industrial areas.

This kit is available in 5 case options: a Robust rugged plastic case, Versatile Softpack reinforced fabric bag, a Weather-resistant plastic



RESPONDER 4 SERIES PLASTIC WATERPROOF FIRST AID KIT 36 X 8.5 X 25CM

AFAK4W

case and a Neat plastic case all our cases and bags are available with custom labelling.

This First Aid Kit comes fully stocked and as part of our <u>Surefill™ instant refill program</u> you can refill and service your client's kits quickly with refill pouches designed to fit into each case type. They provide quick access to all product through a unique centre slit pocket without the need of unpacking each individual product.

CONTENTS LISTS

Code	Description	Qty
AP10125	AEROPLAST Plastic Standard Strip 7.2 x 1.9cm Env/25	25
AFP50125	AEROPLAST Premium Fabric Standard Strip 7.2 x 1.9cm Env/25	25
AF50	AEROFORM Conforming Bandage 5cm x 4M	3
AF75	AEROFORM Conforming Bandage 7.5cm x 4M	3
AFH100	AEROFORM Heavyweight Conforming Bandage 10cm x 4M	1
ABC10	AEROBAND Cotton Calico Triangular Bandage 110 x 110 x 155cm	2
APD50S	AEROPAD Non-Adherent Dressing 5 x 5cm	6
APD100S	AEROPAD Non-Adherent Dressing 7.5 x 10cm	3
APD101S	AEROPAD Non-Adherent Dressing 10 x 10cm	2
AEP1S	AEROPAD Non-Adherent Eye Pads 5.5 x 7.7cm	4
ACD1020S	AEROWOUND Combine Dressing 10 x 20cm	1
AWD14S	AEROWOUND #14 Wound Dressing 12 x 12cm	1
AWD15S	AEROWOUND #15 Wound Dressing 18 x 18cm	1
AGS753S	AEROSWAB Sterile White Gauze Swab 7.5 x 7.5cm (Packs of 3)	5
AP325	AEROTAPE White Microporous Paper Tape 2.5cm x 5M	1
AB3	AEROBURN Burn Gel Sachet 3.5g	5
AAS50	AEROAID First Aid Spray 50ml	1
AW1000	AEROWASH Sodium Chloride Eyewash Ampoule 15ml	8
AW8000	AEROWIPE Alcohol-Free Cleansing Wipes	10
AFS001	AEROSHIELD Disposable Face Shield with one-way valve	1
AGNPF02	AEROGLOVE Large Nitrile Powder-Free Gloves	5
AII2500	AEROCOOL Instant Ice Pack 80g	1
ATB130	AERORESCUE Thermal Blanket 140 x 210cm	1
ASN096P	AEROSUPPLIES First Aid Notebook with Pen	1
AAB100	AEROSUPPLIES Amputated Parts Bags Env/3	3
ASC13	AEROINSTRUMENT Stainless Steel Sharp/Blunt Scissors 13cm	1
AF13	AEROINSTRUMENT Stainless Steel Sharp Forceps 13cm	1
ASP100	AEROPROBE Splinter Probes 3.7cm	10
ASP12	AEROSUPPLIES Medium Safety Pins Bag/12	12
CWD001	AEROGUIDE CPR Flow Chart Card 10 x 14cm	1
AGL001	AEROGUIDE First Aid Leaflet	1