



ADDITIONAL INFORMATION

Weight	0.0833 kg
Dimensions	5.5 × 8.5 × 16 cm
Carton Qty	60
Pricing Unit	Unit/1
Size	0.5L
Min Order Qty	1
ARTG ID	N/A
GMDN Code	N/A
UNSPSC Code	42192802
Tax Schedule	Taxable
Is Medical Device	No
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

AEROHAZARD Sharps Disposal Container is designed for dependable containment of used needles, syringes and other clinical sharps in professional, community and home care settings. A rigid polypropylene body with substantive wall thickness offers robust puncture resistance, helping protect users and surfaces from accidental injury. Its construction also supports leak resistance, so fluids remain contained during handling and interim storage between collections.

The lid assembly is engineered for secure closure. Waste is introduced through a guided entry, then the lid locks down to help prevent reopening and discourage access. This straightforward, one-way approach simplifies safe disposal at the point of care while supporting compliance in environments ranging from treatment rooms and vaccination hubs to outreach and mobile services.

The range is compatible with dedicated wire brackets and mobile trolleys, enabling flexible positioning to suit workflow and available space. Clear identification facilitates quick selection in multi-product environments, and the durable container structure withstands routine transport and handling demands. Built around the core requirements of strength, simplicity and value, it delivers reliable day-to-day performance for sharps management.



AEROHAZARD SHARPS DISPOSAL CONTAINER

500ML

SD500

- **Puncture-resistant build:** rigid polypropylene walls help resist penetration from sharps.
- **Leak-resistant design:** secure closure helps contain fluids during handling.
- **Locking lid:** discourages access after disposal for safer interim storage.
- **Workflow ready:** compatible with wire brackets and trolleys for flexible placement.

The **Size 500 mL** option balances capacity and portability, providing enough volume for moderate sharps usage while remaining easy to position on benches, trolleys or in treatment areas where space is limited.

This sharps container is supplied as a single 500 mL unit, giving a compact yet practical capacity for moderate sharps volumes. It suits benchtops, trolleys and vehicle kits where space is limited but dependable puncture and leak-resistant containment with a locking lid is still essential.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Sharps Disposal Container