



ADDITIONAL INFORMATION

Weight	0.76 kg
Dimensions	25 × 25 × 24 cm
Carton Qty	10
Pricing Unit	Unit/1
Size	12.5L
Min Order Qty	1
ARTG ID	N/A
GMDN Code	N/A
UNSPSC Code	42192802
Tax Schedule	Taxable
Is Medical Device	No
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

AEROHAZARD Sharps Disposal Container is designed for dependable containment of used needles, syringes and other clinical sharps in professional, community and home care settings. A rigid polypropylene body with substantive wall thickness offers robust puncture resistance, helping protect users and surfaces from accidental injury. Its construction also supports leak resistance, so fluids remain contained during handling and interim storage between collections.

The lid assembly is engineered for secure closure. Waste is introduced through a guided entry, then the lid locks down to help prevent reopening and discourage access. This straightforward, one-way approach simplifies safe disposal at the point of care while supporting compliance in environments ranging from treatment rooms and vaccination hubs to outreach and mobile services.

The range is compatible with dedicated wire brackets and mobile trolleys, enabling flexible positioning to suit workflow and available space. Clear identification facilitates quick selection in multi-product environments, and the durable container structure withstands routine transport and handling demands. Built around the core requirements of strength, simplicity and value, it delivers reliable day-to-day performance for sharps management.

- **Puncture-resistant build:** rigid polypropylene walls help resist penetration from sharps.
- **Leak-resistant design:** secure closure helps contain fluids during handling.
- **Locking lid:** discourages access after disposal for safer interim storage.
- **Workflow ready:** compatible with wire brackets and trolleys for flexible placement.

The **Size: 12.5 L** option delivers generous capacity for higher-volume sharps disposal in treatment rooms, vaccination hubs and



AEROHAZARD SHARPS DISPOSAL CONTAINER

12.5L

SD12500

workplace clinics, reducing the frequency of container changeover while retaining a manageable footprint.

The 12.5 litre capacity provides ample room for higher-volume sharps disposal while still fitting comfortably on standard benchtops and mobile trolleys. This size helps reduce how often containers need to be changed in busy clinics, vaccination hubs and workplace health centres, supporting smoother workflows.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Sharps Disposal Container