



### ADDITIONAL INFORMATION

<b>Weight</b>	0.94 kg
<b>Dimensions</b>	25 × 25 × 34 cm
<b>Carton Qty</b>	10
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	17.5L
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	N/A
<b>GMDN Code</b>	N/A
<b>UNSPSC Code</b>	42192802
<b>Tax Schedule</b>	Taxable
<b>Is Medical Device</b>	No
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

### PRODUCT DESCRIPTION

AEROHAZARD Sharps Disposal Container is designed for dependable containment of used needles, syringes and other clinical sharps in professional, community and home care settings. A rigid polypropylene body with substantive wall thickness offers robust puncture resistance, helping protect users and surfaces from accidental injury. Its construction also supports leak resistance, so fluids remain contained during handling and interim storage between collections.

The lid assembly is engineered for secure closure. Waste is introduced through a guided entry, then the lid locks down to help prevent reopening and discourage access. This straightforward, one-way approach simplifies safe disposal at the point of care while supporting compliance in environments ranging from treatment rooms and vaccination hubs to outreach and mobile services.

The range is compatible with dedicated wire brackets and mobile trolleys, enabling flexible positioning to suit workflow and available space. Clear identification facilitates quick selection in multi-product environments, and the durable container structure withstands routine transport and handling demands. Built around the core requirements of strength, simplicity and value, it delivers reliable day-to-day performance for sharps management.

- **Puncture-resistant build:** rigid polypropylene walls help resist penetration from sharps.
- **Leak-resistant design:** secure closure helps contain fluids during handling.
- **Locking lid:** discourages access after disposal for safer interim storage.
- **Workflow ready:** compatible with wire brackets and trolleys for flexible placement.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the



# AEROHAZARD SHARPS DISPOSAL CONTAINER

## 17.5L

SD17500

directions for use.

Sharps Disposal Container