

## **AEROHAZARD SHARPS DISPOSAL CONTAINER 2L**

**SD2000** 



## **ADDITIONAL INFORMATION**

Weight 0.176 kg

**Dimensions**  $13 \times 13 \times 19 \text{ cm}$ 

Carton Qty 48

Pricing Unit Unit/1

Size 2L

Min Order Qty 1

ARTG ID N/A
GMDN Code N/A

**UNSPSC** 42192802

Tax Schedule Taxable

Is Medical Device No

**Is Serialized Item** No

Is Dangerous Good No

## **PRODUCT DESCRIPTION**

Aero offers a wide range of waste disposal containers designed to meet user demands for strength, simplicity and value. This range offers outstanding protection against penetrating sharps and leaking fluids due to its rigid polypropyolene construction and wall thickness.

Sharps Disposal Containers are a durable polypropylene sharps disposal container for sharps and infectious or contaminated medical waste. Easy to use, simply insert waste and lock the lid down for secure containment. Available in a variety of sizes to suit your needs. Durable wire brackets and trolleys also available to suit units.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Sharps Disposal Container



## **AEROHAZARD SHARPS DISPOSAL CONTAINER 2L**

SD2000