



ADDITIONAL INFORMATION

Weight	3.2 kg
Dimensions	14.5 × 14 × 27.5 cm
Carton Qty	5
Pricing Unit	Unit/1
Size	2L
Min Order Qty	1
ARTG ID	N/A
GMDN Code	N/A
UNSPSC	42192802
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

Includes 2 x [SD2000](#) inside steel case.

Aero offers a wide range of waste disposal containers designed to meet user demands for strength, simplicity and value. This range offers outstanding protection against penetrating sharps and leaking fluids due to its rigid polypropylene construction and wall thickness.

Sharps Disposal Containers are a durable polypropylene sharps disposal container for sharps and infectious or contaminated medical waste. Easy to use, simply insert waste and lock the lid down for secure containment. Available in a variety of sizes to suit your needs. Durable wire brackets and trolleys also available to suit units.

Steel Sharps Disposal Container



AEROHAZARD STEEL SHARPS DISPOSAL SAFE 2L (INCLUDES 2 X SD2000)

SD2000L