



ADDITIONAL INFORMATION

Weight	8 kg
Dimensions	187 × 6 × 46 cm
Carton Qty	1
Pricing Unit	Unit/1
Size	1840mm x 450mm x 60mm
Min Order Qty	1
ARTG ID	236801
GMDN Code	13818
Medical Device Class	Class 1
UNSPSC Code	42171613
Tax Schedule	Taxable
Item Length	184
Item Width	45
Item Height	6
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

STR-07

Weight limit: 135kg

Spine board Stretchers facilitate transfer and extrication of patients from vehicles and confined spaces.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AERORESCUE PLASTIC SPINE BOARD STRETCHER WITH STRAPS

STR-07

Plastic Spine Board Stretcher