


AERORESCUE PLASTIC SPINE BOARD STRETCHER WITH STRAPS

STR-07



ADDITIONAL INFORMATION

Weight 8 kg

Dimensions $187 \times 6 \times 46 \text{ cm}$

Carton Qty 1

Pricing Unit Unit/1

Size 1840mm x 450mm x 60mm

Min Order Qty 1

UNSPSC 42171613
Tax Schedule Taxable

Item Length184Item Width45Item Height6Item Measurement Unitcm

Is Serialized Item No
Is Dangerous Good No

PRODUCT DESCRIPTION

AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

STR-07

Weight limit: 135kg

Spine board Stretchers facilitate transfer and extrication of patients from vehicles and confined spaces.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



AERORESCUE PLASTIC SPINE BOARD STRETCHER WITH STRAPS

STR-07

Plastic Spine Board Stretcher