



## ADDITIONAL INFORMATION

<b>Weight</b>	16 kg
<b>Carton Qty</b>	1
<b>Pricing Unit</b>	Unit/1
<b>Size</b>	500mmW x 500mmD x 970mmH
<b>Min Order Qty</b>	1
<b>UNSPSC Code</b>	24101643
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	50
<b>Item Width</b>	50
<b>Item Height</b>	97
<b>Item Measurement Unit</b>	cm
<b>Is Serialized Item</b>	No
<b>Is Dangerous Good</b>	No

## PRODUCT DESCRIPTION

### Medical/Instrument Trolleys

These units provide a sturdy surface and storage for medical supplies of instruments. Stainless steel construction makes for a durable, aesthetic and easy to disinfect finish. Supplied flat packed with minor assembly required.

### Product Features:

- 1 Drawer: 125mm deep
- High quality 304 grade stainless steel
- Safety fold to all edges
- Rubber wheels, front 2 locking
- Full extension telescopic ball bearing drawer runners



# SMALL STAINLESS STEEL TROLLEY WITH DRAWER 50 X 50 X 97CM

ACR12

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Trolley with Drawer