



ADDITIONAL INFORMATION

Weight	7 kg
Carton Qty	1
Pricing Unit	Unit/1
Size	500mmW x 500mmD x 970mmH
Min Order Qty	1
UNSPSC Code	24101643
Tax Schedule	Taxable
Item Length	50
Item Width	50
Item Height	97
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

Medical/Instrument Trolleys

These units provide a sturdy surface and storage for medical supplies of instruments. Stainless steel construction makes for a durable, aesthetic and easy to disinfect finish. Supplied flat packed with minor assembly required.

Product Features:

- Guard rails to 3 sides of 2 shelves
- High quality 304 grade stainless steel
- Safety fold to all edges
- Rubber wheels, front 2 locking



SMALL STAINLESS STEEL TROLLEY WITH RAILS 50 X 50 X 97CM

ACR07

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Trolley with Rails