

AERORESCUE ALLOY COLLAPSIBLE STAIR CHAIR

STR-11



ADDITIONAL INFORMATION

Weight 9.26 kg

Dimensions $100 \times 16 \times 54 \text{ cm}$

Carton Qty 1

Pricing Unit Unit/1

Size 1170mm x 610mm x 910mm

Min Order Qty 1

UNSPSC 42192207
Tax Schedule Taxable

Item Length117Item Width61

Item Height 91

Item Measurement Unit cm

Is Serialized Item No

Is Dangerous Good No

PRODUCT DESCRIPTION

AeroRescue Stretchers are a light-weight high strength solution for many and varied patient recovery situations. Available in a variety of designs and materials for use in a wide range of environments and applications.

STR-11

Weight limit: 160kg

Lightweight aluminium alloy construction specially designed for transporting patients up or down stairs. Handles extend for greater leverage.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Alloy Staircase Ambulance Stretcher



AERORESCUE ALLOY COLLAPSIBLE STAIR CHAIR

STR-11