



### ADDITIONAL INFORMATION

<b>Weight</b>	2.64 kg
<b>Dimensions</b>	36 × 13 × 25 cm
<b>Carton Qty</b>	4
<b>Pricing Unit</b>	Kit
<b>Size</b>	36cm W x 14cm D x 23cm H
<b>Min Order Qty</b>	1
<b>ARTG ID</b>	406543
<b>Medical Device Class</b>	Class 2a
<b>Tax Schedule</b>	Taxable
<b>Item Length</b>	36
<b>Item Width</b>	14
<b>Item Height</b>	23
<b>Item Measurement Unit</b>	cm
<b>Location</b>	AU-Cameron Dr
<b>Is Serialized Item</b>	No
<b>Kits - People (High Risk)</b>	10
<b>Kits - People (Low Risk)</b>	25
<b>Kits - Treatments</b>	Burns, CPR, Dressings and Bandages, Eye Injury, Major Wounds, Minor Wounds, Snake Bite
<b>Kits - Industry</b>	Construction / Outdoor / Remote Location, Education – Excursions / Offsite Activities, Education – Grounds Keeping, Education – Office / Classroom / Hall, Education – Sports / Playground, Home / Personal, Office / Retail, Vehicle
<b>Kits - Case</b>	Soft Pack
<b>Is Dangerous Good</b>	No



### PRODUCT DESCRIPTION

Built for workplaces that need speed, clarity and compliance, the Modulator 4 Series Softpack First Aid Kit brings six colour-coded, injury-specific modules together in a durable, easy-to-carry case. Each module is clearly labelled so responders can identify the right treatment at a glance, streamlining response and reducing time to care during incidents. The softpack uses tough ripstop fabric and a wide-opening design to make packing and inspections simple between jobs and shifts.

Developed to meet and exceed standard workplace regulation contents across Australia, this kit supports low-risk sites up to 1-50 people and high-risk sites up to 1-25 people when stocked and maintained appropriately. It includes a snake bite bandage with a built-in tension indicator, helping users achieve correct application pressure during pressure-immobilisation technique scenarios. Common injury areas are covered via dedicated modules for Cuts & Grazes, Eye Wounds, Dressings & Bandages, Burns, Snake Bite and Trauma, with key treatment cards included to guide confident use.

The 36 × 14 × 23 cm soft case is compact enough for vehicles and site cabinets yet large enough to hold a comprehensive selection of first aid supplies. Clear internal layout and colour coding support quick audits and straightforward replenishment—individual module replacement keeps downtime low and service tasks predictable.

- **Workplace-ready clarity:** six colour-coded modules make selection fast
- **Compliance support:** contents align with Safework Code of Practice for Workplaces
- **Includes tension-indicator bandage:** guides correct pressure for snake bite care
- **Durable softpack:** ripstop fabric with easy-open access for rapid response
- **Audit-friendly layout:** clear labelling and module replacement ease servicing
- **Compact footprint:** 36 × 14 × 23 cm case suits vehicles, site rooms and mobile teams

This softpack configuration features a durable, versatile case that houses all six injury specific modules in one organised system, ideal when portability and easy handling are priorities.

Supplied as 1 complete kit containing 6 individual modules ready for placement in workplaces, homes or vehicles.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.