




ADDITIONAL INFORMATION

Weight	1.02 kg
Dimensions	34 × 22.2 × 4.5 cm
Pricing Unit	Each
Min Order Qty	1
UNSPSC	42301501
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

- Replicates gun-shot wound
- Life-cast silicone simulates human skin, adding realism to improve trainee skill and confidence
- Allows for repeated training in wound packing and haemostatic dressing techniques.
- Tubing for bleeding capability, with blood source controlled via a syringe
- Includes syringe and stage blood

TraumaSim Haemostatic Clotting Trainer Laceration



TRAUMASIM GUNSHOT HAEMOSTATIC CLOTTING TRAINER (SQUARE)

HCT04A