

## TRAUMASIM ARM HAEMORRHAGE CONTROL STABILISING BRACKET

HCS01A



#### **ADDITIONAL INFORMATION**

Weight 0.98 kg

**Dimensions**  $24.1 \times 11.2 \times 19.7 \text{ cm}$ 

Pricing Unit Each

Min Order Qty 1

**UNSPSC** 42301501

Tax Schedule Taxable

Is Serialized Item No
Is Dangerous Good No



#### **PRODUCT DESCRIPTION**

This bracket system can be used to secure the TraumaSim Haemorrhage Control Lifecast Arm (HCS01) arm to a bench or table. 2x additional drill holes are available should you wish to permanently secure the bracket and/or arm to a sturdy structure – additional bolts & washers are required for this, not supplied.

Includes:

Sturdy & durable plastic bracket

2 x 30mm hex-head bolts

2 x cast iron g-clamps

### Check out the HCS01TraumaSim Haemorrhage Control Simulator Arm product here

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

TraumaSim Haemorrhage Control Simulator Arm



# TRAUMASIM ARM HAEMORRHAGE CONTROL STABILISING BRACKET

HCS01A