



ADDITIONAL INFORMATION

Weight	0.98 kg
Dimensions	24.1 × 11.2 × 19.7 cm
Pricing Unit	Each
Min Order Qty	1
UNSPSC	42301501
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

This bracket system can be used to secure the TraumaSim Haemorrhage Control Lifecast Arm (HCS01) arm to a bench or table. 2x additional drill holes are available should you wish to permanently secure the bracket and/or arm to a sturdy structure – additional bolts & washers are required for this, not supplied.

Includes:

Sturdy & durable plastic bracket

2 x 30mm hex-head bolts

2 x cast iron g-clamps

[Check out the HCS01 TraumaSim Haemorrhage Control Simulator Arm product here](#)



TRAUMASIM ARM HAEMORRHAGE CONTROL STABILISING BRACKET

HCS01A