



ADDITIONAL INFORMATION

Weight	11 kg
Dimensions	64 × 49 × 36 cm
Pricing Unit	Each
Min Order Qty	1
Tax Schedule	Taxable
Is Medical Device	No
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The Haemorrhage Control Simulator is a realistic, robust training aid for teaching and practising the most critical and difficult of external haemorrhages – junctional wounds. This training aid with 3 separate wounds – groin wound (junctional), outer thigh wound & mid-thigh amputation – allows trainees to use a junctional tourniquet, a standard limb tourniquet and haemostatic dressings all in one unit. The fleshy wound end can also be used for the stump dressing technique. The simulated pelvis connects to a pump, and simulated blood flows from the wounds; taps control the 3 injuries separately. Set consists of: Includes: Simulated pelvis & thigh All tubing and fittings Silicone skin & urethane bone 5-litre manual pressure pump Rugged construction for field use Packed in rigid plastic trunk Please note connection colours may vary depending on stock availability.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.