

TRAUMASIM HAEMOSTATIC CLOTTING TRAINER GUNSHOT WOUND

HCT03A



ADDITIONAL INFORMATION

Weight 1.18 kg

Dimensions $30 \times 22 \times 5 \text{ cm}$

Pricing Unit Each

Min Order Qty 1

Tax Schedule Taxable

Is Medical Device No
Is Serialized Item No
Is Dangerous Good No

PRODUCT DESCRIPTION

The Haemostatic Clotting Trainer range allows users to pack realistic bleeding wounds to practice the vital skill of haemostatic dressing application. Life-cast silicone simulates human skin, adding realism to improve trainee skill and confidence. Our task trainers allow students to train repeatedly in wound packing and haemostatic dressing techniques with gauze via a sturdy, realistic silicone base with tubing for bleeding capability. The rounded units are to better reflect the curve of the human body (arm or thigh) and are also available as a Gunshot Wound or Laceration with the blood source controlled via a bottle. This unit is a critical training aid for any professional delivering lifesaving first aid training.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

TraumaSim Haemostatic Clotting Trainer Laceration