



ADDITIONAL INFORMATION

Weight	1.18 kg
Dimensions	30 × 22 × 5 cm
Pricing Unit	Each
Min Order Qty	1
Tax Schedule	Taxable
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

The Haemostatic Clotting Trainer range allows users to pack realistic bleeding wounds to practice the vital skill of haemostatic dressing application. Life-cast silicone simulates human skin, adding realism to improve trainee skill and confidence. Our task trainers allow students to train repeatedly in wound packing and haemostatic dressing techniques with gauze via a sturdy, realistic silicone base with tubing for bleeding capability. The rounded units are to better reflect the curve of the human body (arm or thigh) and are also available as a Gunshot Wound or Laceration with the blood source controlled via a bottle. This unit is a critical training aid for any professional delivering lifesaving first aid training.

TraumaSim Haemostatic Clotting Trainer Laceration