



ADDITIONAL INFORMATION

Weight	1.006 kg
Dimensions	33 × 84 × 8.5 cm
Pricing Unit	Each
Min Order Qty	1
UNSPSC Code	42301501
Tax Schedule	Taxable
Is Medical Device	No
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

This medically accurate part-task trainer is made from medical-grade silicone. Improve trainee skills and confidence under pressure using the Emergency Tourniquet Trainer. The original (and world's first!) tourniquet trainer.

The Haemostatic Clotting Trainer range allows users to pack realistic bleeding wounds to practice the vital skill of haemostatic dressing application. Life-cast silicone skin looks and feels real adding realism to improve trainee skill and confidence. These units have the ability to simulate blood flowing from the wound using a pump bottle via tubing. The blood flowing from the wound is stopped when the wound is correctly dressed.

The rounded units are designed to better reflect the curve of the human body (arm or thigh) and Laceration with the blood source controlled via a bottle.

Australian designed and manufactured, this simulator is a necessity for any Trainers looking for suitable practical training aids – and is guaranteed to add realism to your course delivery!

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.



TRAUMASIM LACERATION HAEMOSTATIC CLOTTING TRAINER

HCT03