

ADDITIONAL INFORMATION

Weight	0.04 kg
Dimensions	15 × 0.9 × 6.3 cm
Carton Qty	180
Pricing Unit	Scissor
Size	15cm
Min Order Qty	15
ARTG ID	236810
GMDN Code	13481
Medical Device Class	Class 1
UNSPSC Code	42291615
Tax Schedule	Taxable
Item Length	15
Item Measurement Unit	cm
Is Medical Device	Yes
Is Serialized Item	No
Is Dangerous Good	No



PRODUCT DESCRIPTION

AeroInstrument shears in the AEROINSTRUMENTS range deliver dependable cutting performance for first aid applications. Designed within the Instruments : Scissors category and identified as Trauma Shears, this stainless steel build provides consistent strength for routine kit use when rapid removal of garments or materials is required. The device is comprised of two pivoted blades with finger and thumb holes for controlled handling and is intended for cutting materials such as clothing during initial care. It is not intended for surgical use.

Stainless Steel Construction Manufactured from high-quality stainless steel to support durability, edge integrity and easy cleaning between uses as permitted by local protocols. The metal finish assists smooth cutting while retaining a professional appearance in

clinical, workplace and field kits.

Functional Design The 19 cm form factor offers practical reach and leverage. Finger and thumb rings provide a secure, familiar grip for left- or right-hand techniques. Assorted sizes and specifications are available across the AEROINSTRUMENTS portfolio to match common first aid procedures, helping standardise kit contents for training and workplace compliance.

Compliance & Labelling Products in this line are presented with key handling and regulatory symbols, including CE and UKCA marking, Keep Dry and Keep Away from Sunlight handling cues, plus Non-Sterile and Medical Device identifiers, aiding quick recognition and correct storage. Country of origin details and sponsor/manufacturer addresses are provided on pack for transparent supply chain documentation.

Pack & Identification Clear code and description panels (e.g., USP19 – Stainless Steel Universal Shears, 19 cm / 7.5 in) support accurate reordering and inventory checks. Barcode labelling and consistent AeroInstrument branding streamline scanning and stock control for resellers and distributors. The overall presentation is optimised for inclusion in first aid kits, emergency response bags and workplace cabinets where reliable, easy-to-identify instruments are essential.

The 15 cm configuration with a plastic protective tip offers compact control for safe garment removal in first aid. Stainless steel blades provide reliable cutting action, and the kit-friendly size fits neatly in workplace, vehicle and home kits.

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

Universal Scissors