



ADDITIONAL INFORMATION

Weight	0.26 kg
Dimensions	18 × 15 × 8 cm
Carton Qty	25
Pricing Unit	Unit/1
Min Order Qty	1
Tax Schedule	Taxable
Is Serialized Item	No
Kits - People (Low Risk)	5
Kits - Treatments	Burns, CPR, Eye Injury, Major Wounds, Minor Wounds
Kits - Industry	Home, Vehicle
Kits - Case	Refill
Is Dangerous Good	No
Barcode (GTIN/EAN/UPC)	9341394010116



PRODUCT DESCRIPTION

CONTENTS LISTS

Code	Description	Qty
AP1011	AEROPLAST Plastic Plasters 72 mm x 19 mm	10
AF50	AEROFORM Conforming Bandage 5 cm x 4 m	1
AF75	AEROFORM Conforming Bandage 7.5 cm x 4 m	1
ABC10	AEROBAND Calico Triangular Bandage 110 cm x 110 cm x 155 cm	1
APD101S	AEROPAD Non-Adherent Dressing 10 cm x 10 cm	2
AWD14S	AEROWOUND Wound Dressing #14 Medium BPC	1
AGS753S	AEROSWAB Gauze Swabs 7.5 cm x 7.5 cm 3's	3
AP325	AEROTAPE Microporous Paper Tape 2.5 cm x 5 m	1
AB3	AEROBURN Burn Gel Sachet 3.5 g	1
AW1000	AEROWASH Eye Wash Ampoule 15 mL	1
AW8000	AEROWIPE Cleansing Wipe Alcohol Free	3
AFS001	AEROSHIELD CPR Face Shield Disposable	1

AGNPF02	AEROGLOVE Nitrile Examination Gloves Pair	1
ASN096P	AEROSUPPLIES Notebook & Pen	1
ASC13	AEROINSTRUMENT Scissors 13 cm	1
AF13	AEROINSTRUMENT Forceps 11 cm	1
ASP12	AEROSUPPLIES Safety Pins Medium	3
AGL001	AEROGUIDE First Aid Leaflet	1