



ADDITIONAL INFORMATION

Weight	0.005 kg
Dimensions	14 × 12 × 4 cm
Carton Qty	100
Pricing Unit	Unit/1
Size	75mm
Min Order Qty	1
ARTG ID	257365
GMDN Code	42424
Medical Device Class	Class 1s
UNSPSC	42271913
Tax Schedule	Non-Taxable
Item Length	7.5
Item Measurement Unit	cm
Is Serialized Item	No
Is Dangerous Good	No

PRODUCT DESCRIPTION

Guedel Airways are a critical item for any advanced trauma response to keep airways clear of obstruction when providing advanced first response treatment.

Guedel Airways are available in a range of sizes which must be appropriately matched to the patient to ensure a correct fit. Colours help with ease of recognition of a suitable size. Each airway is individually packaged for best hygiene practice.

Guedel Airway #0 is colour black, and is 60mm in diameter.

White Guedel Airway Size 1.5



WHITE SMALL STERILE GUEDEL AIRWAY SIZE 1.5 (ISO 7.0)

GAB1015