
# HEARTSTART HS1 CHILD SMART PADS PAIR/2

**HS1-PAD-PED**




### ADDITIONAL INFORMATION

| | |
|---|---|
| **Weight** | 0.36 kg |
| **Dimensions** | 21.7 × 4 × 13.9 cm |
| **Pricing Unit** | Each |
| **Tax Schedule** | Non-Taxable |
| **Is Serialized Item** | No |
| **Is Dangerous Good** | No |

### PRODUCT DESCRIPTION

- Child / Infant pads for SCA victims under 8 y.o., or weighing less than 25kg.
- Packaged in a cartridge that can clips into the AED
- Adhesive pads
- Clearly printed diagrams to show rescuers the correct pad placement, as well as sensors to indicate when the pads have been place correctly, improving the effectiveness of the defibrilator's shock.
- 2 year shelf life
- Compatible with M5066A

This product is not a substitute for direct medical attention from a healthcare practitioner. Always read the label and follow the directions for use.

# HEARTSTART HS1 CHILD SMART PADS PAIR/2

**HS1-PAD-PED**